



designed for scientists



C 6040 CalWin

/// Data Sheet

The PC software CalWin C 6040 allows to save, evaluate and easy transfer of even larger volumes of data to your laboratory information management system (LIMS).

Multiple devices can be controlled by the software.

It functions ideally in conjunction with the new C 6000 global standards, the C 6000 isoperibol, C 1, the C 200 (Firmware 2.0) as well as C 5000 (Firmware 2.22) and the C 2000 (Firmware 2.22).

> Clear display of all measurements, results and connected calorimeters



designed for scientists

- > Printing and storage of calibration logs with control chart view
- > Library and grouping function with enhanced data filter options
- > Correction calculations to determine the lower heating value (bet calorific value) according to ISO, ASTM, DIN, GB, GOST and other standards
- > Data export into Excel worksheets or in csv (ASCII code) files

Scope of delivery:

- > DVD with software and PDF instructions for installation
- > 0003036000 C 5041.10 PC - calorimeter RS 232 connecting cable (3 m)
- > 0004589500 RS 232 - USB adapter

Minimum PC requirements:

- > Windows XP (SP2), Windows Vista, Windows 7, Windows 8, Windows 10
- > at least one free USB connection or one RS 232 (9-pin sub-D (M) serial interface
- > processor min. 1.6 GHz (single core processor)
- > 2 GB RAM
- > 2.5 GB hard disk space
- > DVD-ROM drive



designed for scientists