

AMBINOVA TECHNOLOGIES PVT. LTD.

ISO : 9001:2015

ISO : 13485:2016

VORTEX MIXER



FEATURES

- Press mode operation
- High speed Mixing
- Speed Up to 3500 rpm
- Orbit 4.5 mm
- Continuous Mode Operation
- With cup attachment
- DC motor.
- ABS plastic body.

MODEL NO.	DESCRIPTION
AMBI-VM-01	Ambinova Vortex Mixer cup attachment
AMBI-VM-03	Foam attachment for 1.5 and 2 ml tubes

ABOUT COMPANY

Ambinova Technologies Pvt Ltd is dynamic and forward-thinking manufacturing company specializing in laboratory instruments designed for Life Science and Healthcare Laboratories. With a commitment to excellence, Ambinova has earned a reputation for producing high-quality, laboratory instruments that meet the evolving needs of scientific research, diagnostics, and analytical testing.

Mission:

Through innovative design, rigorous quality control, and a commitment to customer satisfaction, we strive to be the premier provider of reliable, accurate, and cutting-edge laboratory instruments. With a dedication to continuous improvement and technological advancement, we aim to enhance the efficiency, productivity, and success of scientific community.

Vision:

Our vision is to be one of the most trusted and value-adding companies Globally.

Products and Services:

Our extensive range of laboratory equipment includes Drybath, Incubators, Autoclaves, Laminar Air Flow, Shakers, and Laboratory Furniture. Each product is meticulously crafted using the qualified and tested materials and cutting-edge technology, ensuring optimal performance and reliability. Our talented team work collaboratively to drive progress, solve complex challenges, and deliver exceptional value to our clients.

CERTIFICATIONS



+91 8595173503
+91 9625922608



sales@ambinova.in



Plot at KH No.9/25/2, 13/5/1 Industrial Area Nangli Sakrawati,
Najafgarh, New Delhi- 110043.